

Accufuser EasyBolus

Model	Basal rate (ml/hr)	Bolus(ml)	Lock-out time (min)	Volume (ml)
E20210M	2	2	10	100
E20215M	2	2	15	100
E20230M	2	2	30	100
E20260M	2	2	60	100
E20310M	3	2	10	100
E20315M	3	2	15	100
E20330M	3	2	30	100
E20360M	3	2	60	100
E20410M	4	2	10	100
E20415M	4	2	15	100
E20430M	4	2	30	100
E20460M	4	2	60	100
E20510M	5	2	10	100
E20515M	5	2	15	100
E20530M	5	2	30	100
E20560M	5	2	60	100
E20810M	8	2	10	100
E20815M	8	2	15	100
E20830M	8	2	30	100
E20860M	8	2	60	100
E21010M	10	2	10	100
E21015M	10	2	15	100
E21030M	10	2	30	100
E21060M	10	2	60	100
E20210L	2	2	10	300
E20215L	2	2	15	300
E20230L	2	2	30	300
E20260L	2	2	60	300
E20310L	3	2	10	300
E20315L	3	2	15	300
E20330L	3	2	30	300
E20360L	3	2	60	300
E20410L	4	2	10	300
E20415L	4	2	15	300

E20430L	4	2	30	300
E20460L	4	2	60	300
E20510L	5	2	10	300
E20515L	5	2	15	300
E20530L	5	2	30	300
E20560L	5	2	60	300
E20810L	8	2	10	300
E20815L	8	2	15	300
E20830L	8	2	30	300
E20860L	8	2	60	300
E21010L	10	2	10	300
E21015L	10	2	15	300
E21030L	10	2	30	300
E21060L	10	2	60	300
E20210XL	2	2	10	600
E20215XL	2	2	15	600
E20230XL	2	2	30	600
E20260XL	2	2	60	600
E20310XL	3	2	10	600
E20315XL	3	2	15	600
E20330XL	3	2	30	600
E20360XL	3	2	60	600
E20410XL	4	2	10	600
E20415XL	4	2	15	600
E20430XL	4	2	30	600
E20460XL	4	2	60	600
E20510XL	5	2	10	600
E20515XL	5	2	15	600
E20530XL	5	2	30	600
E20560XL	5	2	60	600
E20810XL	8	2	10	600
E20815XL	8	2	15	600
E20830XL	8	2	30	600
E20860XL	8	2	60	600
E21010XL	10	2	10	600
E21015XL	10	2	15	600
E21030XL	10	2	30	600

E21060XL	10	2	60	600
----------	----	---	----	-----

Model	Basal rate (ml/hr)	Bolus(ml)	Lock-out time (min)	Volume (ml)
E30210M	2	3	10	100
E30215M	2	3	15	100
E30230M	2	3	30	100
E30260M	2	3	60	100
E30310M	3	3	10	100
E30315M	3	3	15	100
E30330M	3	3	30	100
E30360M	3	3	60	100
E30410M	4	3	10	100
E30415M	4	3	15	100
E30430M	4	3	30	100
E30460M	4	3	60	100
E30510M	5	3	10	100
E30515M	5	3	15	100
E30530M	5	3	30	100
E30560M	5	3	60	100
E30810M	8	3	10	100
E30815M	8	3	15	100
E30830M	8	3	30	100
E30860M	8	3	60	100
E31010M	10	3	10	100
E31015M	10	3	15	100
E31030M	10	3	30	100
E31060M	10	3	60	100
E30210L	2	3	10	300
E30215L	2	3	15	300
E30230L	2	3	30	300
E30260L	2	3	60	300
E30310L	3	3	10	300
E30315L	3	3	15	300
E30330L	3	3	30	300
E30360L	3	3	60	300
E30410L	4	3	10	300
E30415L	4	3	15	300

E30430L	4	3	30	300
E30460L	4	3	60	300
E30510L	5	3	10	300
E30515L	5	3	15	300
E30530L	5	3	30	300
E30560L	5	3	60	300
E30810L	8	3	10	300
E30815L	8	3	15	300
E30830L	8	3	30	300
E30860L	8	3	60	300
E31010L	10	3	10	300
E31015L	10	3	15	300
E31030L	10	3	30	300
E31060L	10	3	60	300
E30210XL	2	3	10	600
E30215XL	2	3	15	600
E30230XL	2	3	30	600
E30260XL	2	3	60	600
E30310XL	3	3	10	600
E30315XL	3	3	15	600
E30330XL	3	3	30	600
E30360XL	3	3	60	600
E30410XL	4	3	10	600
E30415XL	4	3	15	600
E30430XL	4	3	30	600
E30460XL	4	3	60	600
E30510XL	5	3	10	600
E30515XL	5	3	15	600
E30530XL	5	3	30	600
E30560XL	5	3	60	600
E30810XL	8	3	10	600
E30815XL	8	3	15	600
E30830XL	8	3	30	600
E30860XL	8	3	60	600
E31010XL	10	3	10	600
E31015XL	10	3	15	600
E31030XL	10	3	30	600

E31060XL	10	3	60	600
----------	----	---	----	-----

Model	Basal rate (ml/hr)	Bolus(ml)	Lock-out time (min)	Volume (ml)
E40210M	2	4	10	100
E40215M	2	4	15	100
E40230M	2	4	30	100
E40260M	2	4	60	100
E40310M	3	4	10	100
E40315M	3	4	15	100
E40330M	3	4	30	100
E40360M	3	4	60	100
E40410M	4	4	10	100
E40415M	4	4	15	100
E40430M	4	4	30	100
E40460M	4	4	60	100
E40510M	5	4	10	100
E40515M	5	4	15	100
E40530M	5	4	30	100
E40560M	5	4	60	100
E40810M	8	4	10	100
E40815M	8	4	15	100
E40830M	8	4	30	100
E40860M	8	4	60	100
E41010M	10	4	10	100
E41015M	10	4	15	100
E41030M	10	4	30	100
E41060M	10	4	60	100
E40210L	2	4	10	300
E40215L	2	4	15	300
E40230L	2	4	30	300
E40260L	2	4	60	300
E40310L	3	4	10	300
E40315L	3	4	15	300
E40330L	3	4	30	300
E40360L	3	4	60	300
E40410L	4	4	10	300
E40415L	4	4	15	300

E40430L	4	4	30	300
E40460L	4	4	60	300
E40510L	5	4	10	300
E40515L	5	4	15	300
E40530L	5	4	30	300
E40560L	5	4	60	300
E40810L	8	4	10	300
E40815L	8	4	15	300
E40830L	8	4	30	300
E40860L	8	4	60	300
E41010L	10	4	10	300
E41015L	10	4	15	300
E41030L	10	4	30	300
E41060L	10	4	60	300
E40210XL	2	4	10	600
E40215XL	2	4	15	600
E40230XL	2	4	30	600
E40260XL	2	4	60	600
E40310XL	3	4	10	600
E40315XL	3	4	15	600
E40330XL	3	4	30	600
E40360XL	3	4	60	600
E40410XL	4	4	10	600
E40415XL	4	4	15	600
E40430XL	4	4	30	600
E40460XL	4	4	60	600
E40510XL	5	4	10	600
E40515XL	5	4	15	600
E40530XL	5	4	30	600
E40560XL	5	4	60	600
E40810XL	8	4	10	600
E40815XL	8	4	15	600
E40830XL	8	4	30	600
E40860XL	8	4	60	600
E41010XL	10	4	10	600
E41015XL	10	4	15	600
E41030XL	10	4	30	600

E41060XL	10	4	60	600
----------	----	---	----	-----

Model	Basal rate (ml/hr)	Bolus(ml)	Lock-out time (min)	Volume (ml)
E50210M	2	5	10	100
E50215M	2	5	15	100
E50230M	2	5	30	100
E50260M	2	5	60	100
E50310M	3	5	10	100
E50315M	3	5	15	100
E50330M	3	5	30	100
E50360M	3	5	60	100
E50410M	4	5	10	100
E50415M	4	5	15	100
E50430M	4	5	30	100
E50460M	4	5	60	100
E50510M	5	5	10	100
E50515M	5	5	15	100
E50530M	5	5	30	100
E50560M	5	5	60	100
E50810M	8	5	10	100
E50815M	8	5	15	100
E50830M	8	5	30	100
E50860M	8	5	60	100
E51010M	10	5	10	100
E51015M	10	5	15	100
E51030M	10	5	30	100
E51060M	10	5	60	100
E50210L	2	5	10	300
E50215L	2	5	15	300
E50230L	2	5	30	300
E50260L	2	5	60	300
E50310L	3	5	10	300
E50315L	3	5	15	300
E50330L	3	5	30	300
E50360L	3	5	60	300
E50410L	4	5	10	300
E50415L	4	5	15	300

E50430L	4	5	30	300
E50460L	4	5	60	300
E50510L	5	5	10	300
E50515L	5	5	15	300
E50530L	5	5	30	300
E50560L	5	5	60	300
E50810L	8	5	10	300
E50815L	8	5	15	300
E50830L	8	5	30	300
E50860L	8	5	60	300
E51010L	10	5	10	300
E51015L	10	5	15	300
E51030L	10	5	30	300
E51060L	10	5	60	300
E50210XL	2	5	10	600
E50215XL	2	5	15	600
E50230XL	2	5	30	600
E50260XL	2	5	60	600
E50310XL	3	5	10	600
E50315XL	3	5	15	600
E50330XL	3	5	30	600
E50360XL	3	5	60	600
E50410XL	4	5	10	600
E50415XL	4	5	15	600
E50430XL	4	5	30	600
E50460XL	4	5	60	600
E50510XL	5	5	10	600
E50515XL	5	5	15	600
E50530XL	5	5	30	600
E50560XL	5	5	60	600
E50810XL	8	5	10	600
E50815XL	8	5	15	600
E50830XL	8	5	30	600
E50860XL	8	5	60	600
E51010XL	10	5	10	600
E51015XL	10	5	15	600
E51030XL	10	5	30	600

E51060XL	10	5	60	600
----------	----	---	----	-----