

Accufuser Omnibus

Model	A Mode Flow rate (ml/hr)	B Mode Flow rate (ml/hr)	C Mode Flow rate (ml/hr)	Bolus volume (min)	Lockout Time (min)	Volume (ml)
O2020310M	2	3	5	2	10	100
O2020315M	2	3	5	2	15	100
O2020330M	2	3	5	2	30	100
O2020360M	2	3	5	2	60	100
O2020410M	2	4	6	2	10	100
O2020415M	2	4	6	2	15	100
O2020430M	2	4	6	2	30	100
O2020460M	2	4	6	2	60	100
O2030510M	3	5	8	2	10	100
O2030515M	3	5	8	2	15	100
O2030530M	3	5	8	2	30	100
O2030560M	3	5	8	2	60	100
O2030610M	3	6	9	2	10	100
O2030615M	3	6	9	2	15	100
O2030630M	3	6	9	2	30	100
O2030660M	3	6	9	2	60	100
O2040610M	4	6	10	2	10	100
O2040615M	4	6	10	2	15	100
O2040630M	4	6	10	2	30	100
O2040660M	4	6	10	2	60	100
O2040710M	4	7	11	2	10	100
O2040715M	4	7	11	2	15	100
O2040730M	4	7	11	2	30	100
O2040760M	4	7	11	2	60	100
O2040810M	4	8	12	2	10	100
O2040815M	4	8	12	2	15	100
O2040830M	4	8	12	2	30	100
O2040860M	4	8	12	2	60	100
O2050710M	5	7	12	2	10	100
O2050715M	5	7	12	2	15	100
O2050730M	5	7	12	2	30	100
O2050760M	5	7	12	2	60	100
O2050810M	5	8	13	2	10	100
O2050815M	5	8	13	2	15	100

O2050830M	5	8	13	2	30	100
O2050860M	5	8	13	2	60	100
O2051010M	5	10	15	2	10	100
O2051015M	5	10	15	2	15	100
O2051030M	5	10	15	2	30	100
O2051060M	5	10	15	2	60	100
O2081010M	8	10	18	2	10	100
O2081015M	8	10	18	2	15	100
O2081030M	8	10	18	2	30	100
O2081060M	8	10	18	2	60	100
O2020310L	2	3	5	2	10	300
O2020315L	2	3	5	2	15	300
O2020330L	2	3	5	2	30	300
O2020360L	2	3	5	2	60	300
O2020410L	2	4	6	2	10	300
O2020415L	2	4	6	2	15	300
O2020430L	2	4	6	2	30	300
O2020460L	2	4	6	2	60	300
O2030510L	3	5	8	2	10	300
O2030515L	3	5	8	2	15	300
O2030530L	3	5	8	2	30	300
O2030560L	3	5	8	2	60	300
O2030610L	3	6	9	2	10	300
O2030615L	3	6	9	2	15	300
O2030630L	3	6	9	2	30	300
O2030660L	3	6	9	2	60	300
O2040610L	4	6	10	2	10	300
O2040615L	4	6	10	2	15	300
O2040630L	4	6	10	2	30	300
O2040660L	4	6	10	2	60	300
O2040710L	4	7	11	2	10	300
O2040715L	4	7	11	2	15	300
O2040730L	4	7	11	2	30	300
O2040760L	4	7	11	2	60	300
O2040810L	4	8	12	2	10	300
O2040815L	4	8	12	2	15	300
O2040830L	4	8	12	2	30	300

O2040860L	4	8	12	2	60	300
O2050710L	5	7	12	2	10	300
O2050715L	5	7	12	2	15	300
O2050730L	5	7	12	2	30	300
O2050760L	5	7	12	2	60	300
O2050810L	5	8	13	2	10	300
O2050815L	5	8	13	2	15	300
O2050830L	5	8	13	2	30	300
O2050860L	5	8	13	2	60	300
O2051010L	5	10	15	2	10	300
O2051015L	5	10	15	2	15	300
O2051030L	5	10	15	2	30	300
O2051060L	5	10	15	2	60	300
O2081010L	8	10	18	2	10	300
O2081015L	8	10	18	2	15	300
O2081030L	8	10	18	2	30	300
O2081060L	8	10	18	2	60	300
O2020310XL	2	3	5	2	10	600
O2020315XL	2	3	5	2	15	600
O2020330XL	2	3	5	2	30	600
O2020360XL	2	3	5	2	60	600
O2020410XL	2	4	6	2	10	600
O2020415XL	2	4	6	2	15	600
O2020430XL	2	4	6	2	30	600
O2020460XL	2	4	6	2	60	600
O2030510XL	3	5	8	2	10	600
O2030515XL	3	5	8	2	15	600
O2030530XL	3	5	8	2	30	600
O2030560XL	3	5	8	2	60	600
O2030610XL	3	6	9	2	10	600
O2030615XL	3	6	9	2	15	600
O2030630XL	3	6	9	2	30	600
O2030660XL	3	6	9	2	60	600
O2040610XL	4	6	10	2	10	600
O2040615XL	4	6	10	2	15	600
O2040630XL	4	6	10	2	30	600
O2040660XL	4	6	10	2	60	600

O2040710XL	4	7	11	2	10	600
O2040715XL	4	7	11	2	15	600
O2040730XL	4	7	11	2	30	600
O2040760XL	4	7	11	2	60	600
O2040810XL	4	8	12	2	10	600
O2040815XL	4	8	12	2	15	600
O2040830XL	4	8	12	2	30	600
O2040860XL	4	8	12	2	60	600
O2050710XL	5	7	12	2	10	600
O2050715XL	5	7	12	2	15	600
O2050730XL	5	7	12	2	30	600
O2050760XL	5	7	12	2	60	600
O2050810XL	5	8	13	2	10	600
O2050815XL	5	8	13	2	15	600
O2050830XL	5	8	13	2	30	600
O2050860XL	5	8	13	2	60	600
O2051010XL	5	10	15	2	10	600
O2051015XL	5	10	15	2	15	600
O2051030XL	5	10	15	2	30	600
O2051060XL	5	10	15	2	60	600
O2081010XL	8	10	18	2	10	600
O2081015XL	8	10	18	2	15	600
O2081030XL	8	10	18	2	30	600
O2081060XL	8	10	18	2	60	600
O3020310M	2	3	5	3	10	100
O3020315M	2	3	5	3	15	100
O3020330M	2	3	5	3	30	100
O3020360M	2	3	5	3	60	100
O3020410M	2	4	6	3	10	100
O3020415M	2	4	6	3	15	100
O3020430M	2	4	6	3	30	100
O3020460M	2	4	6	3	60	100
O3030510M	3	5	8	3	10	100
O3030515M	3	5	8	3	15	100
O3030530M	3	5	8	3	30	100
O3030560M	3	5	8	3	60	100
O3030610M	3	6	9	3	10	100

O3030615M	3	6	9	3	15	100
O3030630M	3	6	9	3	30	100
O3030660M	3	6	9	3	60	100
O3040610M	4	6	10	3	10	100
O3040615M	4	6	10	3	15	100
O3040630M	4	6	10	3	30	100
O3040660M	4	6	10	3	60	100
O3040710M	4	7	11	3	10	100
O3040715M	4	7	11	3	15	100
O3040730M	4	7	11	3	30	100
O3040760M	4	7	11	3	60	100
O3040810M	4	8	12	3	10	100
O3040815M	4	8	12	3	15	100
O3040830M	4	8	12	3	30	100
O3040860M	4	8	12	3	60	100
O3050710M	5	7	12	3	10	100
O3050715M	5	7	12	3	15	100
O3050730M	5	7	12	3	30	100
O3050760M	5	7	12	3	60	100
O3050810M	5	8	13	3	10	100
O3050815M	5	8	13	3	15	100
O3050830M	5	8	13	3	30	100
O3050860M	5	8	13	3	60	100
O3051010M	5	10	15	3	10	100
O3051015M	5	10	15	3	15	100
O3051030M	5	10	15	3	30	100
O3051060M	5	10	15	3	60	100
O3081010M	8	10	18	3	10	100
O3081015M	8	10	18	3	15	100
O3081030M	8	10	18	3	30	100
O3081060M	8	10	18	3	60	100
O3020310L	2	3	5	3	10	300
O3020315L	2	3	5	3	15	300
O3020330L	2	3	5	3	30	300
O3020360L	2	3	5	3	60	300
O3020410L	2	4	6	3	10	300
O3020415L	2	4	6	3	15	300

O3020430L	2	4	6	3	30	300
O3020460L	2	4	6	3	60	300
O3030510L	3	5	8	3	10	300
O3030515L	3	5	8	3	15	300
O3030530L	3	5	8	3	30	300
O3030560L	3	5	8	3	60	300
O3030610L	3	6	9	3	10	300
O3030615L	3	6	9	3	15	300
O3030630L	3	6	9	3	30	300
O3030660L	3	6	9	3	60	300
O3040610L	4	6	10	3	10	300
O3040615L	4	6	10	3	15	300
O3040630L	4	6	10	3	30	300
O3040660L	4	6	10	3	60	300
O3040710L	4	7	11	3	10	300
O3040715L	4	7	11	3	15	300
O3040730L	4	7	11	3	30	300
O3040760L	4	7	11	3	60	300
O3040810L	4	8	12	3	10	300
O3040815L	4	8	12	3	15	300
O3040830L	4	8	12	3	30	300
O3040860L	4	8	12	3	60	300
O3050710L	5	7	12	3	10	300
O3050715L	5	7	12	3	15	300
O3050730L	5	7	12	3	30	300
O3050760L	5	7	12	3	60	300
O3050810L	5	8	13	3	10	300
O3050815L	5	8	13	3	15	300
O3050830L	5	8	13	3	30	300
O3050860L	5	8	13	3	60	300
O3051010L	5	10	15	3	10	300
O3051015L	5	10	15	3	15	300
O3051030L	5	10	15	3	30	300
O3051060L	5	10	15	3	60	300
O3081010L	8	10	18	3	10	300
O3081015L	8	10	18	3	15	300
O3081030L	8	10	18	3	30	300

O3081060L	8	10	18	3	60	300
O3020310XL	2	3	5	3	10	600
O3020315XL	2	3	5	3	15	600
O3020330XL	2	3	5	3	30	600
O3020360XL	2	3	5	3	60	600
O3020410XL	2	4	6	3	10	600
O3020415XL	2	4	6	3	15	600
O3020430XL	2	4	6	3	30	600
O3020460XL	2	4	6	3	60	600
O3030510XL	3	5	8	3	10	600
O3030515XL	3	5	8	3	15	600
O3030530XL	3	5	8	3	30	600
O3030560XL	3	5	8	3	60	600
O3030610XL	3	6	9	3	10	600
O3030615XL	3	6	9	3	15	600
O3030630XL	3	6	9	3	30	600
O3030660XL	3	6	9	3	60	600
O3040610XL	4	6	10	3	10	600
O3040615XL	4	6	10	3	15	600
O3040630XL	4	6	10	3	30	600
O3040660XL	4	6	10	3	60	600
O3040710XL	4	7	11	3	10	600
O3040715XL	4	7	11	3	15	600
O3040730XL	4	7	11	3	30	600
O3040760XL	4	7	11	3	60	600
O3040810XL	4	8	12	3	10	600
O3040815XL	4	8	12	3	15	600
O3040830XL	4	8	12	3	30	600
O3040860XL	4	8	12	3	60	600
O3050710XL	5	7	12	3	10	600
O3050715XL	5	7	12	3	15	600
O3050730XL	5	7	12	3	30	600
O3050760XL	5	7	12	3	60	600
O3050810XL	5	8	13	3	10	600
O3050815XL	5	8	13	3	15	600
O3050830XL	5	8	13	3	30	600
O3050860XL	5	8	13	3	60	600

O3051010XL	5	10	15	3	10	600
O3051015XL	5	10	15	3	15	600
O3051030XL	5	10	15	3	30	600
O3051060XL	5	10	15	3	60	600
O3081010XL	8	10	18	3	10	600
O3081015XL	8	10	18	3	15	600
O3081030XL	8	10	18	3	30	600
O3081060XL	8	10	18	3	60	600
O4020310M	2	3	5	4	10	100
O4020315M	2	3	5	4	15	100
O4020330M	2	3	5	4	30	100
O4020360M	2	3	5	4	60	100
O4020410M	2	4	6	4	10	100
O4020415M	2	4	6	4	15	100
O4020430M	2	4	6	4	30	100
O4020460M	2	4	6	4	60	100
O4030510M	3	5	8	4	10	100
O4030515M	3	5	8	4	15	100
O4030530M	3	5	8	4	30	100
O4030560M	3	5	8	4	60	100
O4030610M	3	6	9	4	10	100
O4030615M	3	6	9	4	15	100
O4030630M	3	6	9	4	30	100
O4030660M	3	6	9	4	60	100
O4040610M	4	6	10	4	10	100
O4040615M	4	6	10	4	15	100
O4040630M	4	6	10	4	30	100
O4040660M	4	6	10	4	60	100
O4040710M	4	7	11	4	10	100
O4040715M	4	7	11	4	15	100
O4040730M	4	7	11	4	30	100
O4040760M	4	7	11	4	60	100
O4040810M	4	8	12	4	10	100
O4040815M	4	8	12	4	15	100
O4040830M	4	8	12	4	30	100
O4040860M	4	8	12	4	60	100
O4050710M	5	7	12	4	10	100

O4050715M	5	7	12	4	15	100
O4050730M	5	7	12	4	30	100
O4050760M	5	7	12	4	60	100
O4050810M	5	8	13	4	10	100
O4050815M	5	8	13	4	15	100
O4050830M	5	8	13	4	30	100
O4050860M	5	8	13	4	60	100
O4051010M	5	10	15	4	10	100
O4051015M	5	10	15	4	15	100
O4051030M	5	10	15	4	30	100
O4051060M	5	10	15	4	60	100
O4081010M	8	10	18	4	10	100
O4081015M	8	10	18	4	15	100
O4081030M	8	10	18	4	30	100
O4081060M	8	10	18	4	60	100
O4020310L	2	3	5	4	10	300
O4020315L	2	3	5	4	15	300
O4020330L	2	3	5	4	30	300
O4020360L	2	3	5	4	60	300
O4020410L	2	4	6	4	10	300
O4020415L	2	4	6	4	15	300
O4020430L	2	4	6	4	30	300
O4020460L	2	4	6	4	60	300
O4030510L	3	5	8	4	10	300
O4030515L	3	5	8	4	15	300
O4030530L	3	5	8	4	30	300
O4030560L	3	5	8	4	60	300
O4030610L	3	6	9	4	10	300
O4030615L	3	6	9	4	15	300
O4030630L	3	6	9	4	30	300
O4030660L	3	6	9	4	60	300
O4040610L	4	6	10	4	10	300
O4040615L	4	6	10	4	15	300
O4040630L	4	6	10	4	30	300
O4040660L	4	6	10	4	60	300
O4040710L	4	7	11	4	10	300
O4040715L	4	7	11	4	15	300

O4040730L	4	7	11	4	30	300
O4040760L	4	7	11	4	60	300
O4040810L	4	8	12	4	10	300
O4040815L	4	8	12	4	15	300
O4040830L	4	8	12	4	30	300
O4040860L	4	8	12	4	60	300
O4050710L	5	7	12	4	10	300
O4050715L	5	7	12	4	15	300
O4050730L	5	7	12	4	30	300
O4050760L	5	7	12	4	60	300
O4050810L	5	8	13	4	10	300
O4050815L	5	8	13	4	15	300
O4050830L	5	8	13	4	30	300
O4050860L	5	8	13	4	60	300
O4051010L	5	10	15	4	10	300
O4051015L	5	10	15	4	15	300
O4051030L	5	10	15	4	30	300
O4051060L	5	10	15	4	60	300
O4081010L	8	10	18	4	10	300
O4081015L	8	10	18	4	15	300
O4081030L	8	10	18	4	30	300
O4081060L	8	10	18	4	60	300
O4020310XL	2	3	5	4	10	600
O4020315XL	2	3	5	4	15	600
O4020330XL	2	3	5	4	30	600
O4020360XL	2	3	5	4	60	600
O4020410XL	2	4	6	4	10	600
O4020415XL	2	4	6	4	15	600
O4020430XL	2	4	6	4	30	600
O4020460XL	2	4	6	4	60	600
O4030510XL	3	5	8	4	10	600
O4030515XL	3	5	8	4	15	600
O4030530XL	3	5	8	4	30	600
O4030560XL	3	5	8	4	60	600
O4030610XL	3	6	9	4	10	600
O4030615XL	3	6	9	4	15	600
O4030630XL	3	6	9	4	30	600

O4030660XL	3	6	9	4	60	600
O4040610XL	4	6	10	4	10	600
O4040615XL	4	6	10	4	15	600
O4040630XL	4	6	10	4	30	600
O4040660XL	4	6	10	4	60	600
O4040710XL	4	7	11	4	10	600
O4040715XL	4	7	11	4	15	600
O4040730XL	4	7	11	4	30	600
O4040760XL	4	7	11	4	60	600
O4040810XL	4	8	12	4	10	600
O4040815XL	4	8	12	4	15	600
O4040830XL	4	8	12	4	30	600
O4040860XL	4	8	12	4	60	600
O4050710XL	5	7	12	4	10	600
O4050715XL	5	7	12	4	15	600
O4050730XL	5	7	12	4	30	600
O4050760XL	5	7	12	4	60	600
O4050810XL	5	8	13	4	10	600
O4050815XL	5	8	13	4	15	600
O4050830XL	5	8	13	4	30	600
O4050860XL	5	8	13	4	60	600
O4051010XL	5	10	15	4	10	600
O4051015XL	5	10	15	4	15	600
O4051030XL	5	10	15	4	30	600
O4051060XL	5	10	15	4	60	600
O4081010XL	8	10	18	4	10	600
O4081015XL	8	10	18	4	15	600
O4081030XL	8	10	18	4	30	600
O4081060XL	8	10	18	4	60	600
O5020310M	2	3	5	5	10	100
O5020315M	2	3	5	5	15	100
O5020330M	2	3	5	5	30	100
O5020360M	2	3	5	5	60	100
O5020410M	2	4	6	5	10	100
O5020415M	2	4	6	5	15	100
O5020430M	2	4	6	5	30	100
O5020460M	2	4	6	5	60	100

O5030510M	3	5	8	5	10	100
O5030515M	3	5	8	5	15	100
O5030530M	3	5	8	5	30	100
O5030560M	3	5	8	5	60	100
O5030610M	3	6	9	5	10	100
O5030615M	3	6	9	5	15	100
O5030630M	3	6	9	5	30	100
O5030660M	3	6	9	5	60	100
O5040610M	4	6	10	5	10	100
O5040615M	4	6	10	5	15	100
O5040630M	4	6	10	5	30	100
O5040660M	4	6	10	5	60	100
O5040710M	4	7	11	5	10	100
O5040715M	4	7	11	5	15	100
O5040730M	4	7	11	5	30	100
O5040760M	4	7	11	5	60	100
O5040810M	4	8	12	5	10	100
O5040815M	4	8	12	5	15	100
O5040830M	4	8	12	5	30	100
O5040860M	4	8	12	5	60	100
O5050710M	5	7	12	5	10	100
O5050715M	5	7	12	5	15	100
O5050730M	5	7	12	5	30	100
O5050760M	5	7	12	5	60	100
O5050810M	5	8	13	5	10	100
O5050815M	5	8	13	5	15	100
O5050830M	5	8	13	5	30	100
O5050860M	5	8	13	5	60	100
O5051010M	5	10	15	5	10	100
O5051015M	5	10	15	5	15	100
O5051030M	5	10	15	5	30	100
O5051060M	5	10	15	5	60	100
O5081010M	8	10	18	5	10	100
O5081015M	8	10	18	5	15	100
O5081030M	8	10	18	5	30	100
O5081060M	8	10	18	5	60	100
O5020310L	2	3	5	5	10	300

O5020315L	2	3	5	5	15	300
O5020330L	2	3	5	5	30	300
O5020360L	2	3	5	5	60	300
O5020410L	2	4	6	5	10	300
O5020415L	2	4	6	5	15	300
O5020430L	2	4	6	5	30	300
O5020460L	2	4	6	5	60	300
O5030510L	3	5	8	5	10	300
O5030515L	3	5	8	5	15	300
O5030530L	3	5	8	5	30	300
O5030560L	3	5	8	5	60	300
O5030610L	3	6	9	5	10	300
O5030615L	3	6	9	5	15	300
O5030630L	3	6	9	5	30	300
O5030660L	3	6	9	5	60	300
O5040610L	4	6	10	5	10	300
O5040615L	4	6	10	5	15	300
O5040630L	4	6	10	5	30	300
O5040660L	4	6	10	5	60	300
O5040710L	4	7	11	5	10	300
O5040715L	4	7	11	5	15	300
O5040730L	4	7	11	5	30	300
O5040760L	4	7	11	5	60	300
O5040810L	4	8	12	5	10	300
O5040815L	4	8	12	5	15	300
O5040830L	4	8	12	5	30	300
O5040860L	4	8	12	5	60	300
O5050710L	5	7	12	5	10	300
O5050715L	5	7	12	5	15	300
O5050730L	5	7	12	5	30	300
O5050760L	5	7	12	5	60	300
O5050810L	5	8	13	5	10	300
O5050815L	5	8	13	5	15	300
O5050830L	5	8	13	5	30	300
O5050860L	5	8	13	5	60	300
O5051010L	5	10	15	5	10	300
O5051015L	5	10	15	5	15	300

O5051030L	5	10	15	5	30	300
O5051060L	5	10	15	5	60	300
O5081010L	8	10	18	5	10	300
O5081015L	8	10	18	5	15	300
O5081030L	8	10	18	5	30	300
O5081060L	8	10	18	5	60	300
O5020310XL	2	3	5	5	10	600
O5020315XL	2	3	5	5	15	600
O5020330XL	2	3	5	5	30	600
O5020360XL	2	3	5	5	60	600
O5020410XL	2	4	6	5	10	600
O5020415XL	2	4	6	5	15	600
O5020430XL	2	4	6	5	30	600
O5020460XL	2	4	6	5	60	600
O5030510XL	3	5	8	5	10	600
O5030515XL	3	5	8	5	15	600
O5030530XL	3	5	8	5	30	600
O5030560XL	3	5	8	5	60	600
O5030610XL	3	6	9	5	10	600
O5030615XL	3	6	9	5	15	600
O5030630XL	3	6	9	5	30	600
O5030660XL	3	6	9	5	60	600
O5040610XL	4	6	10	5	10	600
O5040615XL	4	6	10	5	15	600
O5040630XL	4	6	10	5	30	600
O5040660XL	4	6	10	5	60	600
O5040710XL	4	7	11	5	10	600
O5040715XL	4	7	11	5	15	600
O5040730XL	4	7	11	5	30	600
O5040760XL	4	7	11	5	60	600
O5040810XL	4	8	12	5	10	600
O5040815XL	4	8	12	5	15	600
O5040830XL	4	8	12	5	30	600
O5040860XL	4	8	12	5	60	600
O5050710XL	5	7	12	5	10	600
O5050715XL	5	7	12	5	15	600
O5050730XL	5	7	12	5	30	600

O5050760XL	5	7	12	5	60	600
O5050810XL	5	8	13	5	10	600
O5050815XL	5	8	13	5	15	600
O5050830XL	5	8	13	5	30	600
O5050860XL	5	8	13	5	60	600
O5051010XL	5	10	15	5	10	600
O5051015XL	5	10	15	5	15	600
O5051030XL	5	10	15	5	30	600
O5051060XL	5	10	15	5	60	600
O5081010XL	8	10	18	5	10	600
O5081015XL	8	10	18	5	15	600
O5081030XL	8	10	18	5	30	600
O5081060XL	8	10	18	5	60	600